### NEXTSEQ 550 Siyu Sun

## References

- NextSeq 550 System:
  - <u>https://support.illumina.com/sequencing/sequencing\_instruments/nextseq-550/documentation.html</u>
  - <u>https://support.illumina.com/content/dam/illumina-support/documents/documentation/</u> system\_documentation/nextseq-550dx/nextseq-550dx-instrument-ref-guide-1000000009513-07.pdf
- NextSeq 500 and 550 System Denature and Dilute Libraries Guide
  - <u>https://support.illumina.com/content/dam/illumina-support/documents/documentation/</u> system\_documentation/nextseq/15048776\_18\_nextseq-500-550-denature-dilute-libraries-guide.pdf
- PhiX control:
  - <u>https://www.illumina.com/content/dam/illumina-support/documents/products/technotes/</u> <u>technote\_phixcontrolv3.pdf</u>

# before running

- pre thaw the reagent cartridge, HT1 buffer (put in 4C fridge overnight)
- before running, place flow cell, reagent cartridge at RT for at least 30min

- you need to also prepare: wash solution (125ml 0.05% Tween20 wash)
  - you could also use the wash solution next to the machine, but make sure there are enough

## pre step: Wash

- if the previous run on the machine is less than 2 weeks, then  $\bullet$ this step is not necessary;
- if the previous run is more than 2 weeks ago, the machine will prompt you to do a wash step:
  - select "quick wash", then follow the instructions on screen
  - prepare 125ml 0.05% Tween20 wash (usually they have a bottle of wash solution on the side, but you could also bring your own)
  - load the wash solution to the buffer wash cartridge
  - place the buffer wash cartridge (with wash solution),  $\bullet$ reagent wash cartridge, and the empty waste collecting cartridge to place
  - click start, it will take ~20min.

wash buffer load here



### step1: create your run on "Local Run Manager" software

- prepare a sample sheet in .csv format (refer to the sample sheet template), save to your flash drive, plug-in and copy-paste to the desktop
- open "local run manager" software
- log in (the sheet with login name and password info are on the wall to the left)
- click: Create Run—GenerateFASTQ



			kton > VO
NEXTSE NB55164		Organize - New folder	
Create R	UT GenerateFAST	Q Desktop Downloads	
Import Sample S	Sheet	Documents SC	G_SS 230915_SampleS heet.csv
Run Name*		3D Objects	
Run Name		Desktop	
		<ul> <li>Documents</li> <li>Downloads</li> </ul>	
Run Settings		Music	
Library Prep	Select	File name: 2309	915_SampleSheet.csv

- ulletsoftware; or
- you could manually type in the information ullet
- once you make sure every info is there, you click "save run" then minimize the software window ullet

you could select "Import Sample Sheet" to directly upload your pre-generated sample sheet.csv to the



## step2: set-up a sequencing run

- 2.1: Select Run
- in the home screen, select: "Experiment" "Sequence"





- select a run from the list

   select "local Run
   Manager", check
   BaseSpace and
   "Proactive, Run
   monitoring and Storage"
- you will need to log in to your personal illumina account at this point
- you will then input the username and password of "Local Run Manager"





551647				
1				5
Setup	Load	Review	Check	Sequence
t run, monitorin	ng and storage mode.			
ns Using				
Manual				
Local Run Ma	anager			
co Soquence H	ub Settings			
BaseSpace S Increase efficie	equence Hub (requires an ency and reduce errors using	internet connection) online run setup with enhance	⊳ ed validation.	
) Proactive	and Run Monitoring 🕜			
Proactive	, Run Monitoring and Stor	age 🕜		-
Back	🗶 Exit			Next ->

(2)			
Load	Review	Check	Sequence
oring, you must log in to B	aseSpace.		
i	llumina		
EII Address			
Email Address			
Password			
Sign In			
Evit			Next ->

NB551647 ssu	n@uchc.edu			
0				
Run Setup	Load	Review	Check	Sequence
Enter your Local Run N	Manager username and pa	ssword, then select Next.		
			N	
			6	
Usernar	ne* admin			
Passwo	rd *			
	-			
🗲 Back	🗙 Exit			

- created run to the NextSeq RUO software.
- make sure everything is looking correct

<b>NB551647</b> ad	dmin				and the second sec
Image: Constraint of the sector of the se	Load list of runs from Local Run I	(3) Revi Manager. Select	ew Next to Load Co	(4) Check nsumables.	
Q Search for Run Nam	e		Created by Run ID Date Module	Admin Use 2002 2023-09-1 GenerateF	er 5 astQWorkflo
230915_GATA3_PF	10		Read Length Custom Primers Output Folder	Read 1 II 47 D:\Output	ndex 1 Ind
✓ ✓ ✓	C X Exit			Edit	Preview Sa

• after log-in to "local run manager", you are able to extract the previously

• you can see all the info about your experiment here, and noble check to



## 2.2 load flow cell

- make sure the flow cell is placing at RT for at least 30min — so that no condensation generated
- \*\* it will take a while for the machine to automatically open the flow cell placing stage
- remove the used flow cell out
- put the new flow cell on the stage as prompt by machine
- select "Load"
- select "Next"



### 2.3 Empty/replace the Spent Reagent Container

- open the buffer compartment door
- took the container out, empty/clean the container, the waste should be disposed into the <u>sequencing reagent bottle</u> next to the sink
- \*there is a spare one on the shelf, if need it, use it



Figure 13 Remove the Spent Reagents Container





### slide the empty Spent Reagent Container into the buffer compartment (bottom) until it stops





# 2.4 load the buffer cartridge

- remove the used ones  $\bullet$ 
  - \*used ones should be rinsed with water and solution can be disposed into the sink; the container itself should be disposed into regular trash can
- slide the new buffer cartridge into the top stage until it stops
- close the buffer compartment door, and select "Next"



# 2.5 load the reagent cartridge

- reagent cartridge need to pre-thaw, and place in RT for at least 30min before loading
- open the reagent compartment door
- remove the used reagent cartridge
  - used reagent cartridge need to be properly clean-up: step1. remove the rubber cover, press down to eject the reservoir; step2. dispose the solution containing formamide to the bottle next to the sink labeled "Toxic" and "formamide", the reservoir itself goes to the red-biohazard trash can; step3. the remaining cartridge need to be rinsed, and dispose into the regular trash can

Figure 17 Removable Position #6



- Protective rubber cover
- Position #6





Reservoir that are ejectable and contains formamide



- load your library into the reagent cartridge
  - poke the foil with a clean tip, then load with pipette at the designated area
- slide the reagent cartridge into the stage until it stops
- close the door
- select "Load", then "Next"





## step3. monitor run progress

- review check start
- at this point, you will once again review the details of your run, if everything looks good, click "Next"
- machine will then take several minutes to automatically check the temperature, systems ..., once done, click "Start"

NB551647	admin				NB551647 admin			Setting valve to position [, 1]	
O-Run Setup	Load	3 Review	Check	6 Sequence	©- Run Setup	Load	Review	4 Check	6 Sequence
Review the run :	settings.				(i) The instrument is performing	a self-check. When it is	s finished, select Start	to begin sequencing.	
un Name	230915_GATA3_PRO	Module	GenerateFastC	QWorkflow					
escription	230915_GATA3_PRO				System Checks	System Ch	iecks		
lead Type	Paired End			📝 Edit	Temperature	Doors Clo	osed	<ul> <li>Instrument Disk S</li> </ul>	pace
	Read 1 Index 47 6	1Index 2 0	Read 2 <b>39</b>		Imaging System	<ul><li>Consuma</li><li>Required</li></ul>	bles Loaded Software	<ul><li>Network Connecti</li><li>Network Disk Spa</li></ul>	on ce
Custom Primer:					Reagent Delivery	_			
Output Folder	D:\Output Purge consumables for th	iis run		-	Cancel				ß
🗲 Back	🗙 Exit			Next 🛧	Cancer	Exit			Start ->

- you will see the machine are starting to process your library...it will give you an expected time of run completion
  - still need to take 3~4 hours to generate the fasta files for you
- time and your contact, so that if anything happened, you get notified by others



• Notice that this time is the time of the sequencing completion; the machine will

I also recommend placing a sticky note on the machine, write the expected running

Initial Prime: 0 of 1 repeat	s	Running from 9115-9116 Apm. Please Contact	
R2: Cycle 0 of 39		Curtin Lab. Thx 22	
30 Ib %	Clu  Cluste	ers Passing Filter	

## final step: clean-up

- after you done with your run, make sure to clean up everything.
- right place.
- finally, wish your run successful!

make sure you dispose the solutions, the reservoir, the cartridges to the